
**carl
marlin**
Solingen
Finest Dental Instruments

Cervi } ect

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EN - Instructions for use

DE - Gebrauchsanweisung	3-5
EN - Instruction for use	6-8
FR - Mode d'emploi	9-11
NL - Gebruiksaanwijzing	12-14
ES - Instrucciones de uso	15-17
HR - Upute za uporabu	18-20
SV - Bruksanvisning	21-23
DA - Brugsanvisning	24-26
SL - Navodila za uporabo	27-29
NO - Bruksanvisning	30-32
CS - Návod k použití	33-35
IT - Istruzioni per l'uso	36-38
SK - Návod na použitie	39-41
RO - Instrucțiuni de utilizare	42-44
EL - Οδηγίες χρήσης	45-47
LT - Naudojimo instrukcija	48-50
HU - Használati utasítás	51-53
BG - Инструкции за употреба	54-56
PT - Instruções de utilização	57-59
LV - Lietošanas instrukcija	60-62

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Intraligamentäre Spritze CerviJect

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Instructions for use

Intraligamental syringe CerviJect

The injection syringe for intraligamental anaesthesia, for use with a 1.8 ml anaesthesia cartridge and needle with metric (**art. no. 1956**) or US thread (**art. no. 1956-IT**), which can be used in intraseptal anaesthesia.


Intended use only by specialised dental personnel.

Intraligamental anaesthesia (ILA) is a special form of infiltration or terminal anaesthesia. ILA differs from terminal anaesthesia in that the injection is made directly into the gingival sulcus of the tooth in question and not parallel to the tooth axis in the fold of the mucosa next to the teeth to be anaesthetised. This means that only the tooth to be treated can be selectively anaesthetised, whereas terminal anaesthesia usually affects 2-3 teeth.

Function and operation

1. Functional parts and structure of the CerviJect

- | | |
|-------------------------------------|---|
| a. Head section | g. Dosing lever |
| b. Viewing window | h. Spring |
| c. Burst protection | i. Reset button for the plunger |
| d. Bayonet catch | j. Plunger |
| e. Cannula attachment (with thread) | k. Injection needle (not included in the scope of delivery) |
| f. Handpiece | |

2. **CerviJect - processing (disassembled state)**

The CerviJect syringe is processed in its disassembled state:

Head section (a), handpiece (f), plunger (j), burst protection (c). The dosing lever (g) on the handpiece (f) is specially screwed on and must not be removed for processing.

3. **Assembly and preparation for injection**

Insert the plunger (j) into the handpiece (f) as shown and return it completely to the starting position; to do this, press the reset button (i), hold the CerviJect syringe vertically and return the plunger (j) into the handpiece (f) as far as it will go. First insert the burst protection (c) into the head section (a) and then the standard anaesthesia cartridge. Connect the head section (a) and handpiece (f) using the bayonet catch (d).

Operate the dosing lever (g) until slight resistance is felt. The plunger (j) then makes contact with the plug of the anaesthesia cartridge and pushes the cartridge head to the end of the head section (a).

Insert the disposable cannula (k) into the needle attachment (e) and screw it on. We recommend the use of cannulas (the cartridge carrier needles) in the sizes 0.3 x 12 mm or 0.3 x 16 mm with a short bevel corresponding to the anaesthesia method.

CAUTION: Always insert the cartridge into the cartridge holder first and then screw on the cannula attachment (e) (needle). If the order is reversed, there is a risk that the needle will not penetrate the cartridge and will bend. As a result, there is no connection between the needle and the anaesthetic and high pressure is built up in the cartridge when the dosing lever (g) is actuated, so that the cartridge can burst.

Operate the dosing lever (g) until anaesthetic emerges from the needle tip. The syringe is now ready for injection.

4. **Descriptions of use for ILA**

a. **Advantages of ILA for the user**

- Anaesthesia is induced after a short time
- Only a very small injection volume is required
- Targeted, selective anaesthesia of individual teeth
- Elimination of the risk of soft tissue haematomas
- Intracanalally effective e.g. for endodontic treatment
- Also effective for highly acidic tissue/inflammation

b. Advantages for the patient

- Fast and gentle analgesia
- No numbness of neighbouring anatomical structures (lip, tongue)
- Only low intake of anaesthetic (high-risk patients)
- Reduced duration of anaesthesia
- As the ILA syringe does not have the typical appearance of a conventional cartridge syringe, the slim and graceful design has a calming effect on the patient. It can therefore generally also be used successfully on children.

5. Care instructions

The CerviJect syringe is sold non-sterile. Before first use, carry out a complete processing according to our processing instructions "QSA 313" (on our website www.carl martin.de under Downloads).

After each use, loosen the head and handpiece of the syringe and remove the anaesthesia cartridge, burst protection and plunger. After cleaning and disinfection, the individual parts must be dried thoroughly. The individual parts must then be placed in a suitable sterilisation medium (sterile packaging, container) and only opened again shortly before the next use. Only assemble the individual parts when they are dry.

It must be ensured that only validated procedures for cleaning, disinfection and sterilisation are used. The validation of the processing in accordance with our processing instructions "QSA 313" by solgiene oHG confirms that the syringes can be with a standard washer-disinfector in accordance with EN ISO 15883.

Cleaning in the ultrasonic cleaner: The maximum sonication time should not exceed 3 minutes at 35 kHz. Otherwise, follow the manufacturer's instructions.

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Carl Martin GmbH