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Ligaturenschneider

ORTHODONTIA

Important note for the use of Carl Martin Orthodontic pliers with tungsten carbide inserts.

The tungsten carbide inserts used by Carl Martin Solingen are of the highest quality. Nevertheless, there is risk of breakage due to improper use. Please cut the wires only in the maximum specified thickness (indicated on the instrument). The cutting performance and service life of our pliers have been tested with all major types of wire and are at the highest level.

Distal end cutter and side cutter

Please cut the wires only in the maximum specified thickness (indicated on the instrument) and with at least **2mm** away from the working end (see Fig.).

Ligature cutter

Precise cutting of ligatures requires a high degree of sharpness from hard metal blades. Constant use combined with the effect of chemicals used in disinfection and the heat in sterilisation may lead to temporary loss of sharpness. As a result, the wires will no longer be cut with precision, but instead could also become crimped. In such cases, we recommend that these instruments be sharpened.

In addition, we also recommend that the cutting of ligatures should NOT be performed with the tip of the ligature cutter, as this can also result in breakages. We recommend, so do other reputable manufacturers, that ligature wires should be cut about **1 mm** back from the tip of the blades (see Fig.).

Below listed the reasons for possible breakage of the tips:

- When severing a wire right at the tip of the blades, the tip slips back from the wire and the blades close together under high pressure. When this happens, the tips may not be adequately able to absorb the compressive force and break as result
- When severing ligature wire in greater thickness, the cutting edge bears the brunt of the very high pressure applied during the severing operation. Similarly in this case, the compressive force cannot be adequately absorbed. This may lead to breakage.

RISK OF BREAKAGE


MIN. 2MM

RISK OF BREAKAGE

MIN. 2MM

RISK OF BREAKAGE

MIN. 1MM

Distal end cutter

Side cutter

Ligature cutter